

Instructions how to perform an energy calibration in the MAR_aB Φ U environment

R. Lutter

February 8, 2007

Abstract

This document describes how to perform an energy calibration by use of MAR_aB Φ U's
MacroBrowser and the ROOT macro `Encal.C`.

1 How to perform an energy calibration

To perform an energy calibration (gamma or particle) call the `MacroBrowser`:

`MacroBrowser`

A menu will then pop up showing several ROOT macros. Choose `Encal.C` from this list.

You'll get a form (fig. 1) to specify parameters needed to start the calibration:

- `Histogram file (.root)`
Click on the folder button and choose the ROOT file containing your calibration spectra. Then select some histograms from the list below either by **range** (first ... last) or one by one (**single**). Press `apply` to accept selection, `clear` to reset list.
- `Calibration source`
Choose `Co60` or `Eu152` for gammas, `TripleAlpha` for particles
- `Calibration energies`
Enter file containing calibration energies. Default is `$MARABOU/data/encal/energies.dat`. See 2.3 for a format description.
- `Calibration data`
Where to write calibration data (extension `.cal`).
Calibration data will be formatted according to ROOT's resource format (ROOT object `TEnv`, 2.1);
- `Calibration results`
Where to write detailed calibration results (extension `.res`). Same formatting rules as for `.cal` files apply (2.2).
- `Fit results`
ROOT file where to store histograms together with fit data as well as linear regression data (extension `.root`).
- `Clear output files`
Click `yes` to remove existing files (`.cal` and `.res`) on start. If `no` is chosen new data will be appended to existing files, existing items will be overwritten with new values.
- `Precalibration file`
To get an `Eu152` calibration you have to preform a `Co60` calibration step first in order to assign gamma energies to peak centers. The name of the `Co60` calibration file from a previous run has to be given here.
- `Lower/Upper limit`
Limits of the display/fit region (channels)
- `Threshold for peak finder`
Relative height of peaks taken into account given by percentage of maximum peak (continued on next page)

- **Sigma**
sigma value for gaussian fit
- **Fit mode**
Either `gauss` or `gauss+tail`: how to fit lines in the histogram.
- **Fit grouping**
You may fit `single` peaks or an `ensemble` of peaks (may give better results for `Co60` and `TripleAlpha` calibrations).
- **Background**
You have the choice between a `const` or `linear` background.
- **Range for peak fit**
Set lower and upper limits for the peak fit region (single peak or group fitting), given in units of `sigma`.
- **Display results**
Check button `step` to stop after each fit, check `2dim` to display a two-dim histogram containing re-calibrated histograms line by line.
- **Verbose**
Choosing button `yes` will output some diagnostics to `stdout`.




Figure 1: Enca1.C: how to perform an energy calibration

Press **Execute** to start the calibration.

A new window will pop up showing fit results (fig. 2).

The upper part shows peak centroids (**blue markers**) as returned by the peak finder together with the gaussian fit (**red curve**).

The lower part shows the linear regression results when assigning alpha lines 5.157 MeV (^{239}Pu), 5.486 MeV (^{241}Am), and 5.865 MeV (^{244}Cm) to peaks, respectively.

The title bar below the graphs shows the resulting calibration formula $E(x) = a_0 + a_1 * x$.




Figure 2: Encal.C: fit and calibration results

If you have chosen **step** mode you now have to press one of the buttons at the bottom edge of the display:

- **accept** add calibration data to output file (.cal) and continue with next histogram.
- **discard** discard calibration and continue with next histogram.
- **stop** leave calibration loop and return to MacroBrowser's menu (fig. 1).
- **quit** close files and exit from ROOT.

If you checked the `2dim` button you will get a two-dimensional histogram showing calibrated histograms line by line (fig. 3).

From this plot you can see immediately that calibration of histogram `hRawAs1` didn't give satisfactory results. Now try the following:

- Go back to `MacroBrowser`'s menu (fig. 1)
- Reset the histogram list by pressing `clear`
- Select histogram `hRawAs1` from the list (button `apply`)
- Change flag `Clear output files` to `no`
- Change the `sigma` value to `2`
- Press `Execute` to start over

The new fit data for histogram `hRawAs1` will then overwrite the existing ones without changing other entries. You then end up with fig. 4.




Figure 3: Encal .C: 2-dim plot of calibrated histograms, first try




Figure 4: Encal .C: 2-dim plot after re-fitting

2 File formats

2.1 Calibration data (extension .cal)

A calibration data file (extension .cal) generated by Encal.C is formatted according to ROOT's resource format (see ROOT documentation, object TEnv). It consists of a header followed by entries for each histogram.

Header	Calib.ROOTFile: <HistoFile>.root Calib.Source: Co60 or Eu152 or TripleAlpha Calib.Energies: <path to file containing calib energies> Calib.NofHistograms: <NH>
Entry (1 per histo)	Calib.<HistoName>.Xmin: <lower limit> Calib.<HistoName>.Xmax: <upper limit> Calib.<HistoName>.Gain: <gain a_1 > Calib.<HistoName>.Offset: <offset a_0 >

This file may be read in directly to the data acquisition program by use of method TMrbAnalyze::ReadCalibrationFromFile(

2.2 Detailed calibration results (extension .res)

In addition to the calibration data file (2.1) detailed fit results will be written to a file with extension .res. Its format is similar to that of the .cal file but contains much more information.

Header	Calib.ROOTFile: <HistoFile>.root Calib.Source: Co60 or Eu152 or TripleAlpha Calib.Energies: <path to file containing calib energies> Calib.NofHistograms: <NH> Calib.Emin: <minimum energy over all histos> Calib.Emax: <maximum energy> Calib.Xmin: <minimum channel> Calib.Xmax: <maximum channel>
Histo entry (1 per histo)	Calib.<HistoName>.Xmin: <lower limit> Calib.<HistoName>.Xmax: <upper limit> Calib.<HistoName>.NofPeaks: <NP> Calib.<HistoName>.Gain: <gain a_1 > Calib.<HistoName>.Offset: <offset a_0 > Calib.<HistoName>.FitOk: TRUE, FALSE, or AUTO
Peak entry (1 per peak per histo)	Calib.<HistoName>.Peak.<I>.X: <centroid peak finder> Calib.<HistoName>.Peak.<I>.Xfit: <centroid gaussian fit> Calib.<HistoName>.Peak.<I>.Xerr: <centroid error> Calib.<HistoName>.Peak.<I>.Y: <amplitude peak finder> Calib.<HistoName>.Peak.<I>.Yfit: <amplitude gaussian fit> Calib.<HistoName>.Peak.<I>.Yerr: <amplitude error> Calib.<HistoName>.Peak.<I>.Chi2: <chisquare gaussian fit>

2.3 Calibration energies

Calibration energies have to be provided as ROOT resources (see ROOT documentation, object TEnv).

Header	Calib.NofCalibrations: <NC> Calib.SourceNames: <Source1>:<Source2>:...:<SourceNC>
Entry (1 per line per source)	Calib.<SourceN>.NofLines: <NL> Calib.<SourceN>.Line.<I>.E: <energy> Calib.<SourceN>.Line.<I>.Eerr: <error>